

# Multi-Functional Combination Instrument

- ⦿ Vacuum Negative Pressure
- ⦿ 6MHz radio frequency
- ⦿ 180° mechanical rotation
- ⦿ LED red and blue light
- ⦿ Real-time Temperature Control

Face  
Lifting

Anti-  
Wrinkle

Fat  
Burning

Body  
Shaping

Tendons  
relaxing



**PHYSIOESTHÉTIKA**  
— TECHNOLOGIE —



## ROTARY NEGATIVE PRESSURE RF MACHINE PRINCIPLE

The Rotary Negative Pressure RF Machine uses **vacuum negative pressure + 180° mechanical rotation + temperature control radio frequency + LED red and blue light technology** to pull up the skin and subcutaneous muscle tissue, massage the body evenly, dissolve fat and decompose fat cells into fatty acids and glycerin, and then which can be excreted through the lymphatic circulation, eliminate local stubborn fat, while stimulating the contraction and regeneration of collagen fibers in the dermis, reducing cellulite, and making the skin firmer and smoother.

- ◎ Blood circulation
- ◎ Lymphatic detoxification function

**400%** ↗






## INNER BALL ROLLER MACHINE PRINCIPLE

non-invasive mechanical compression micro-vibration + infrared treatment

Through the 360° rotation and rolling of the silicone ball along the roller, the compression micro-vibration is generated.


As the ball rotates and exerts pressure on the skin, the effect of "pulsating compression" is generated, and the reciprocating movement of continuous pushing, pulling and kneading are achieved.

The tissue will experience some pressure and lifting actions without squeezing or damaging the skin.




### Reshaping action

Reshape contours and smooths out wrinkles




### Muscle action

Tone muscles and tightens tissues




### Sedation

Improve oxygen and eliminate pain



### Drain action

Eliminate toxins and help lymphatic drainage



### Vascular effect

Improve blood circulation and nourish tissues.





# BODY SWIRL HANDLE

## Fat melting + Body shaping

- ⊙ vacuum negative pressure
- ⊙ 6MHz RF
- ⊙ 180°mechanical rotation
- ⊙ LED red and blue dual light



### Swivel Head - Large

- Suitable for waist, abdomen, thighs, back.



### Swivel Head - Small

- Suitable for arms and calves.



5 Inch  
Touch Screen

## FACE HANDLE

Effectively anti-wrinkle

Firm and lift skin

- ⊙ vacuum negative pressure
- ⊙ radio frequency heating

2.8 Inch  
Touch Screen



# BODY ROLLER HANDLE

**Tone and firm the skin tissue**

**Reduce cellulite and shape contour**


- ⊙ Mechanical compression  
micro-vibration
- ⊙ Infrared ray technology

- ① **Display Screen**
- ② **Display Speed Gear**
- ③ **Adjust Speed Button**
- ④ **Switch Rotation  
Direction Button**




## Silicone Large Head

- Silicone is soft, comfortable and good for lifting




## Metal - Gold Plated Large Head

- Metal is good for lymphatic.



## Metal - Silver Plated Large Head

- Metal is good for lymphatic.




# FACIAL ROLLER HANDLE

**Reducing swelling**

**Resisting tissue sagging**

**Rejuvenating the face**

- ⦿ Mechanical compression  
micro-vibration
- ⦿ Infrared ray technology




## Silicone Small Head

▪ Silicone is soft, comfortable and good for lifting



## Metal - Gold Plated Small Head

▪ Metal is good for lymphatic.



## Metal - Silver Plated Small Head

▪ Metal is good for lymphatic.

## ① Display Speed Gear

# Five core technologies

## ◎ 6MHz RF

The RF energy heats cellulite, removes fat deposits, and stimulates the contraction of collagen fibers for better blood flow to all tissues.

## ◎ Real-time Temperature Control

Preset the target temperature before treatment, always control the temperature at the set temperature, and monitor and display it in real time.

## ◎ Vacuum Negative Pressure

Through the rotation movement and contact with the vacuum negative pressure at the same time, it stimulates the lymph and blood circulation, accelerates the elimination of edema, eliminates lymphatic obstruction, and enhances blood flow to the tissue.

## ◎ 180° Rotation and Tweaking

180° alternating rotation and twisting stimulate the metabolism in cells, reduce edema, expel toxins, relieve muscle spasms and pain syndromes.

## ◎ LED red and blue double light

LED red and blue dual light irradiates the skin, which can increase the temperature of fat, stimulate the activity of deep cells, promote blood circulation and decompose triglyceride more efficiently, and help it to be excreted.

MHz





# EIGHT ADVANTAGES OF THE PRODUCT



## ◎ **Dual Interface Design**

The body handle has a 5-inch color touch screen, and the face handle has a 2.8-inch color touch screen, which is more intuitive and easier to control and adjust parameters.



## ◎ **Patented Liposuction Handle**

Ergonomically design, with excellent grip and convenient operation, it makes the beautician feel more relaxed and less prone to fatigue and strain.



## ◎ **Non-invasive**

There is no recovery period, and you can do whatever you want, without affecting your normal work and life.



## ◎ **Long-term Effects**

Long-lasting, after treatment, the skin will improve every day due to the continuous production of new collagen, and will achieve more significant and satisfactory results in about 4-6 months.



## ◎ **Intelligent Temperature Control System**

Set the target temperature of the skin before treatment, and always keep the temperature at the set temperature, which can effectively avoid tissue damage.



## ◎ **Real-time Temperature Detection System**

Real-time control of output energy, real-time detection of skin temperature, real-time display of skin temperature during treatment, safe and efficient.



## ◎ **Safe & Efficient**

Purely physical, no damage to the skin, the best and deepest effect.



## ◎ **Comfortable Labor-saving Convenient**

The treatment process is more comfortable, just like manual liposuction to deal with cellulite, which greatly reduces the pain of treatment.

# FIVE FUNCTIONS OF THE PRODUCT

## ⊙ Anti- Wrinkle

Stimulating the contraction and re-generation of collagen fibers in the dermis, and making the skin firmer and smoother.



## ⊙ Fat Burning

Dissolve fat and decompose fat cells into fatty acids and glycerin.



## ⊙ Body Shaping

Create beautiful body curves.



## ⊙ Tendon relaxing

Dredge meridians and promote blood circulation.



## ⊙ Face Lifting

Tightens sagging face/body skin




# EFFECT COMPARISON CHART

◎ **FACE** anti-wrinkle/face lifting/face shaping



◎ **ABDOMEN** Fat burning/ body shaping



◎ **NECK** anti-wrinkle/ lifting




◎ **BUTTOCK** Fat burning/ body shaping



◎ **LEG** Fat burning/ body shaping



◎ **ARM** Fat burning/ body shaping








# TECHNICAL PARAMETER

Product Name	S11
Touch screen	15 inch
Radio frequency	6MHz
Vacuum Negative Pressure	10-80KPA
Angle of rotation	10-180°
Photon	Red and blue double light
Inner ball roller handle speed	675rpm
operating voltage	AC110V~220V±10%, 50Hz-60Hz
power	10-400W
Airbox size	63×76×118cm
Gross weight	About 79.3kg



# BIG TOUCH LCD SCREEN

-  Energy Adjustment
-  Pulse Width Adjustment
-  Frequency Adjustment
-  Real-time Temperature Monitoring

 2.8" TOUCH LCD SCREEN

 5" TOUCH LCD SCREEN

 15" TOUCH LCD SCREEN

